



10/12 Lead Patient Cables

45502 15-pin D-sub Connector w/ thumbscrews	Bircher	365, 375 CARDIOTRACER 375-11
	Nihon Kohden	ECG-5151A, 5403A ECG-6151A, 6306A, 6353A ECG-6551A, 6543A, 7514A ECG-8110A, 8240A, 8330A ECG-9340, 8350A, 8370A
E45503 15-pin D-sub connector slide-in-groove molding	Burdick Mfg p/n 007517 or 007514 007705, 007706	EK 10, Elite E350, E350 I, E550, E560, Eclipse LE, LE I
	Nihon Kohden	Cardiofax ECG 6551 ECG 6511 ECG 6951D ECG 6511S ECG 6851C ECG 8110P
45505 12-pin round connector	Physio Control	LifePak 12
45506 w/ resistors 15-pin D-sub connector w/ thumbscrews	Hewlett Packard HP / Agilent Part # M2461A M3702C	Pagewriter (new version) M1770 Pagewriter 200i M1771 Pagewriter 100 M1772 Pagewriter 200 300PI (newest Cardiograph)
45507 w/ 10k resistors 15-pin D-sub conn. w/ thumbscrews Formerly P/N 45502-BR	Marquette/ Hellige	EK-51
	Schiller	CS-6 AT-2 AT-10
455006	GE Medical Systems / Marquette	12-Lead Multilink Cable for GE Eagle, Solar, Dash Monitors, Cath Lab Systems
45508 15-Pin D-sub with bottom flanges	Burdick New Style 007704, 007725 012-0844-00	Eclipse 4, 400, 8, 800, 850, 850+ LE II, Atria 3000

