



02 REGULATORS (CLICK STYLE)

- Available in standard (284MA) and mini (285 MA) versions
- Aluminum body, brass core
- 0-4, 0-8, 0-15 and 0-25 LPM models
- CGA 870 (Yoke) or CGA 540 (Nut & nipple) inlet connection
- Hose barb outlet standard; DISS connections available
- Locking ring on T-Handle



Catalog Number	Inlet	Outlet	Flow Range	Flow Settings
284MA-8L	CGA 540	Hose Barb	0-8 LPM	0, 0.5, 1, 1.5, 2, 2.5, 3, 4, 5, 6, 7, 8
284MA-15L	CGA 540	Hose Barb	0-15 LPM	0, 0.5, 1, 2, 3, 4, 5, 6, 7, 8, 10, 15
284MA-15L1	CGA 540	Hose Barb + 1 DISS	0-15 LPM	0, 0.5, 1, 2, 3, 4, 5, 6, 7, 8, 10, 15
284MA-25L1	CGA 540	Hose Barb + 1 DISS	0-25 LPM	0, 0.5, 1, 2, 3, 4, 5, 6, 8, 10, 15, 25
284MA-25L2	CGA 540	Hose Barb + 2 DISS	0-25 LPM	0, 0.5, 1, 2, 3, 4, 5, 6, 8, 10, 15, 25
284MA-8LY	CGA 870	Hose Barb	0-8 LPM	0, 0.5, 1, 1.5, 2, 2.5, 3, 4, 5, 6, 7, 8
284MA-15LY	CGA 870	Hose Barb	0-15 LPM	0, 0.5, 1, 2, 3, 4, 5, 6, 7, 8, 10, 15
284MA-15LY1	CGA 870	Hose Barb + 1 DISS	0-15 LPM	0, 0.5, 1, 2, 3, 4, 5, 6, 7, 8, 10, 15
284MA-25LY1	CGA 870	Hose Barb + 1 DISS	0-25 LPM	0, 0.5, 1, 2, 3, 4, 5, 6, 8, 10, 15, 25
284MA-25LY2	CGA 870	Hose Barb + 2 DISS	0-25 LPM	0, 0.5, 1, 2, 3, 4, 5, 6, 8, 10, 15, 25
MINI Click-style Oxygen Regulators				
285MA-8LY	CGA 870	Hose Barb	0-8 LPM	0, 0.5, 1, 1.5, 2, 2.5, 3, 4, 5, 6, 7, 8
285MA-15LY	CGA 870	Hose Barb	0-15 LPM	0, 0.5, 1, 2, 3, 4, 5, 6, 7, 8, 10, 15