










ECG CABLES & LEAD WIRES

| 10/12 LEAD PATIENT CABLES | | | |
|--|--|---|---|
| 45502 15-pin D-sub Connector with thumbscrews | Bircher | 365, 375 CARDIOTRACER 375-11 |  |
| | Nihon Kohden | ECG-5151A, 5403A ECG-6151A, 6306A, 6353A ECG-6551A, 6543A, 7514A ECG-8110A, 8240A, 8330A ECG-9340, 8350A, 8370A | |
| E45503 15-pin D-sub Connector slide in- -groove molding | Burdick Mfg p/n 007517 or 007514 007705, 007706 | EK 10, Elite E350, E350 I, E550, E560, Eclipse LE, LE I |  |
| | Nihon Kohden | Cardiofax ECG 6551 ECG 6511 ECG 6951D ECG 6511S ECG 6851C ECG 8110P | |
| 45505 12-pin round connector | Physio Control | LifePak 12 |  |
| 45506 w/resistors 15-pin D-sub connector w/thumbscrews | Hewlett Packard HP / Agilent Part # M2461A M3702C | Pagewriter (new version) M1770 Pagewriter 200i M1771 Pagewriter 100 M1772 Pagewriter 200 300PI (newest Cardiograph) |  |
| 45507 w/10k resistors 15-pin D-sub conn. w/thumbscrews Formerly P/N 45502-BR | Marquette/Hellige | EK-51 |  |
| | Schiller | CS-6 AT-2 AT-10 | |
| 455006 | GE Medical Systems / Marquette | 12-Lead Multilink Cable for GE Eagle, Solar, Dash Monitors, Cath Lab Systems |  |
| 45508 15-pin D-sub with bottom flanges | Burdick New Style 007704,007725 012-0844-00 | Eclipse 4, 400, 8 800, 850, 850+ LE II, Atria 3000 |  |

• Add: -S for Snap -P for Pinch -B for Banana