



WALL SUCTION VACUUM REGULATORS

Tenacore repairs wall suction regulators on all models for:

- Amvex
- Boehringer
- Chemetron/Vacutron
- Gentec
- Ohio
- Ohmeda
- Precision Medical
and more

Most hospitals find that it is more cost-effective to repair their Suction regulators rather than replace them.

- Competitive Pricing
- Flat Rate Pricing Available
- Exchange Program Available for unrepairable Regulators

Our testing procedures include an in-depth Physical Inspection focusing on the condition of the:

- Chassis/Housing
- Fitting/Connectors
- Vacuum Regulator
- Mode Selector Switch
- Vacuum Gauge

We adhere to the OEM specifications applied to Leak Testing and Suction Accuracy. All units are Cleaned, Disinfected and Overhauled including replacements of O-rings and Filters.



- Repair Parts Available - see our Replacement Parts Section