

WALL SUCTION REGULATORS

Tenacore wall suction regulators are compact, easy-to-use and include many unique features. Each unit includes soft Side Grips™ that provide users with a secure grip on the unit to prevent unintentional dropping. The Side Grips™ are offered in a range of colors to aid in identifying departmental ownership.

FEATURES

- Analog gauge with $\pm 5\%$ accuracy
- Intermittent timing factory set at 15 seconds on, 8 seconds off
- Adjustable timing (intermittent)
- Color-coded Side Grips™
- Pressure relief valve
- Large, easy-to-read graphics
- Glow-in-the-dark scale and needle

Tenacore wall suction regulators can be customized with your logo and side grip colors to represent the brand of your organization.

MODELS

Pediatric

800 Series	Continuous Intermittent
810 Series	Intermittent Only
820 Series	Continuous Only

Adult

1500 Series	Continuous Intermittent
1510 Series	Intermittent Only
1520 Series	Continuous Only
3000 Series	Continuous-Full

Surgical

7000 Series	Continuous-Full
-------------	-----------------